

# Percutek

## AAA Stent Graft System



Name of Medical Device: Percutek AAA Stent Graft System  
Name of Manufacturer: Beijing Percutek Therapeutics Inc.


Beijing Percutek Therapeutics Inc.  
Building 14, No. 26, Yongwang West Rd, Bioengineering & Medical Industry Base,  
Daxing District, Beijing, China  
www.percutek.com  
+86 10 6125 1354  
info@percutek.com

**SIMPLE  
BEST**


**Proximal Landing Zone  
≥10mm**

**AND**


**Infrarenal Angle  
≤75°**




Mainbody of AAA (PABF)				
Product Code	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Total Covered Length (mm)	OD of Delivery System (Fr)
PABF2012110	20	12	110	16Fr
PABF2212110	22	12	110	16Fr
PABF2412110	24	12	110	18Fr
PABF2614110	26	14	110	18Fr
PABF2814110	28	14	110	18Fr
PABF3014110	30	14	110	20Fr
PABF3214110	32	14	110	20Fr
PABF3414110	34	14	110	20Fr
PABF3614110	36	14	110	20Fr




Limb of AAA (PAIL)				
Product Code	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Total Covered Length (mm)	OD of Delivery System (Fr)
PAIL1610060	16	10	60	14Fr
PAIL1613060	16	13	60	14Fr
PAIL1616060	16	16	60	14Fr
PAIL1620060	16	20	60	16Fr
PAIL1624060	16	24	60	18Fr
PAIL1610080	16	10	80	14Fr
PAIL1613080	16	13	80	14Fr
PAIL1616080	16	16	80	14Fr
PAIL1620080	16	20	80	16Fr
PAIL1624080	16	24	80	18Fr
PAIL1610095	16	10	95	14Fr
PAIL1613095	16	13	95	14Fr
PAIL1616095	16	16	95	14Fr
PAIL1620095	16	20	95	16Fr
PAIL1624095	16	24	95	18Fr
PAIL1610110	16	10	110	14Fr
PAIL1613110	16	13	110	14Fr
PAIL1616110	16	16	110	14Fr
PAIL1620110	16	20	110	16Fr
PAIL1624110	16	24	110	18Fr
PAIL1610125	16	10	125	14Fr
PAIL1613125	16	13	125	14Fr
PAIL1616125	16	16	125	14Fr
PAIL1620125	16	20	125	16Fr
PAIL1624125	16	24	125	18Fr
PAIL1610140	16	10	140	14Fr
PAIL1613140	16	13	140	14Fr
PAIL1616140	16	16	140	14Fr
PAIL1620140	16	20	140	16Fr
PAIL1624140	16	24	140	18Fr
PAIL1610155	16	10	155	14Fr
PAIL1613155	16	13	155	14Fr
PAIL1616155	16	16	155	16Fr
PAIL1620155	16	20	155	16Fr
PAIL1624155	16	24	155	18Fr
Iliac Extension				
PAIL1010080	10	10	80	14Fr
PAIL1313080	13	13	80	14Fr
PAIL2020080	20	20	80	16Fr
PAIL2424080	24	24	80	18Fr



AUI (PAUI)				
Product Code	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Total Covered Length (mm)	OD of Delivery System (Fr)
PAUI1814110	18	14	110	16Fr
PAUI2014110	20	14	110	16Fr
PAUI2214110	22	14	110	16Fr
PAUI2414110	24	14	110	18Fr
PAUI2614110	26	14	110	18Fr
PAUI2814110	28	14	110	18Fr
PAUI3014110	30	14	110	20Fr
PAUI3214110	32	14	110	20Fr
PAUI3414110	34	14	110	20Fr
PAUI3614110	36	14	110	20Fr



Cuff (PACF45)				
Product Code	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Total Covered Length (mm)	OD of Delivery System (Fr)
PACF2020045	20	20	45	16Fr
PACF2222045	22	22	45	16Fr
PACF2424045	24	24	45	18Fr
PACF2626045	26	26	45	18Fr
PACF2828045	28	28	45	18Fr
PACF3030045	30	30	45	20Fr
PACF3232045	32	32	45	20Fr
PACF3434045	34	34	45	20Fr
PACF3636045	36	36	45	20Fr



Cuff (PACF70)				
Product Code	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Total Covered Length (mm)	OD of Delivery System (Fr)
PACF2020070	20	20	70	16Fr
PACF2222070	22	22	70	16Fr
PACF2424070	24	24	70	18Fr
PACF2626070	26	26	70	18Fr
PACF2828070	28	28	70	18Fr
PACF3030070	30	30	70	20Fr
PACF3232070	32	32	70	20Fr
PACF3434070	34	34	70	20Fr
PACF3636070	36	36	70	20Fr